



ClearView™ **3**
Multifocal IOL

Give your patients a ClearView

- **Refractive MIOL**
- **Full Range of Vision**
- **No Concentric Rings**
- **Over 93%** of patients reported **freedom from vision correction***
- **Minimal reports of visual disturbances** deemed clinically insignificant and consistent with those experienced with Monofocal IOLs*

*Data on file

**LENS
TEC**™ 

ClearView 3 Lens Specifications

Optic Type:	Refractive, Aspheric
Add Power:	+3.0D at the IOL plane +2.4D at the spectacle plane
Optic Size:	5.75mm
Length:	11.00mm
Angulation:	0 Degrees
Construction:	1-piece
Material:	Hydrophilic Acrylic (26% water content)
Haptic Style:	Closed loop with four-point fixation
Abbe Value:	61.3 @ 35 degrees
Index of Refraction:	1.456
Model Number:	SBL-3
Power Range:	+10.0D to +36.0D 0.25D available +15.0D to +25.0D 0.50D available +10.0D to +36.0D

Constants

SRK-T:	118.0
Haigis:	α_0 - 0.537 α_1 - 0.333 α_2 - 0.126
Barrett:	LF - 1.36
Holladay 1:	SF - 1.22
Holladay II/ Hoffer Q:	pACD - 4.97

Outcomes You Can Trust



- 0.25D Power Increments
- Industry-Leading Manufacturing Tolerances (0.11D vs 0.40D industry standard)
- Proven Hydrophilic Acrylic Material with 29+ years of clinical use



www.clearviewlens.com



FM 38446.

PKB83 Rev 3 Surgeon